



适用于各种应用的抗化学腐蚀软管和硬管

- PTFE 硬管，尺寸范围为 1/8" (3.2mm) 到 1/16" (1.6mm) O.D.
- 适用于夹管阀的硅胶管和Bio-Chem (C-Flex®) 两类Bio-Chem软管。
- 切割成一定长度，以50英尺为单位长度（软管）和 20 或者 200m 为单位长度（硬管）
- 各种管件均有配套接头（请见接头系统说明书）。

PTFE硬管

完全抗化学腐蚀溶剂的惰PTFE硬管几乎可以适用于任何化学溶剂和各种腐蚀物质的使用上，即使在升高的温度也可胜任。PTFE硬管可以适用于蒸汽，化学试验。可以和所有的“Omnifit 螺纹帽”接口，阀，平底口的接头配套使用，这些即为 Gripper, Omni-Lok™ 或者到锥形接头。

硅胶软管

Bio-Chem的Silicone Select™ 硅胶管是铂金弯曲加工，是为高纯度的医用设备所制造，符合以下各种标准验证：USP XXII, Class VI, FDA , USDA 标准和列入FDA 精加工名单。这种软管都是为最严格的夹管阀使用而造，并且适用于各类夹管阀使用。

软管 (C-Flex)

Bio-Chem (C-Flex®) 软管适用于医疗级高纯度液体使用，符合USP XXII, Class VI, FDA 以及USDA标准，并列入FDA 精加工名单。适用于各类夹管阀使用，尤其适用于有气体渗漏问题的系统。

切管器

手工操作的切割器可以有整齐的切口，可以替换刀片。

货号信息

硬管PTFE（铁氟龙）管型号

型号	规格	圈长度	备注信息
008T16-030-20	1.6mm (1/16") O.D. x 0.3mm (.012") I.D.	20m	3005
008T16-030-200	1.6mm (1/16") O.D. x 0.3mm (.012") I.D.	200m	3006
008T16-050-20	1.6mm (1/16") O.D. x 0.5mm (.02") I.D.	20m	3003
008T16-050-200	1.6mm (1/16") O.D. x 0.5mm (.02") I.D.	200m	3004
008T16-080-20	1.6mm (1/16") O.D. x 0.8mm (.031") I.D.	20m	3001
008T16-080-200	1.6mm (1/16") O.D. x 0.8mm (.031") I.D.	200m	3002
008T16-100-20	1.6mm (1/16") O.D. x 1.0mm (.039") I.D.	20m	3007
008T16-100-200	1.6mm (1/16") O.D. x 1.0mm (.039") I.D.	200m	3008
008T32-150-10	3.2mm (1/8") O.D. x 1.5mm (.059") I.D.	10m	透明/半透明
008T32-150-100	3.2mm (1/8") O.D. x 1.5mm (.059") I.D.	100m	透明/半透明

Tolerances:

Inside diameter: $\pm 0.08\text{mm}$ (0.00275")

Outside diameter: $\pm 0.08\text{mm}$ (0.00275")

软管硅胶管

型号	规格	圈长度	备注信息
10025-10S	0.093" (2.36mm) O.D. x 0.010" (0.25mm) I.D.	50 feet	10025-010093S
10025-23S	0.093" (2.36mm) O.D. x 0.023" (0.58mm) I.D.	50 feet	10025-01-23S
10025-50S	0.134" (3.40mm) O.D. x 0.050" (1.27mm) I.D.	50 feet	10025-01-50S
10025-01S	3/32" (2.38mm) O.D. x 1/32" (0.79mm) I.D.	50 feet	-
10025-02S	1/8" (3.17mm) O.D. x 1/16" (1.59mm) I.D.	50 feet	-
10025-03S	3/16" (4.76mm) O.D. x 1/16" (1.59mm) I.D.	50 feet	-
10025-04S	1/4" (6.35mm) O.D. x 3/16" (4.76mm) I.D.	50 feet	-
10025-05S	1/4" (6.35mm) O.D. x 1/8" (3.17mm) I.D.	50 feet	-
10025-06S	3/8" (9.52mm) O.D. x 1/4" (6.35mm) I.D.	50 feet	-
10025-07S	7/16" (11.1mm) O.D. x 5/16" (7.94mm) I.D.	50 feet	-

软管 (C-Flex)

型号	规格	圈长度	备注信息
10025-23B	0.093" (2.36mm) O.D. x 0.023" (0.58mm) I.D.	50 feet	10025-01-23B
10025-01B	3/32" (2.38mm) O.D. x 1/32" (0.79mm) I.D.	50 feet	-
10025-02B	1/8" (3.17mm) O.D. x 1/16" (1.59mm) I.D.	50 feet	-
10025-03B	3/16" (4.76mm) O.D. x 1/16" (1.59mm) I.D.	50 feet	-
10025-05B	1/4" (6.35mm) O.D. x 1/8" (3.17mm) I.D.	50 feet	-
10025-06B	3/8" (9.52mm) O.D. x 1/4" (6.35mm) I.D.	50 feet	-
10025-07B	7/16" (11.1mm) O.D. x 5/16" (7.94mm) I.D.	50 feet	-

切管器

型号	规格
3062	切管器
3063	单独替换刀片



Trademarks:

C-Flex® is a registered trademark of Consolidated Polymer Technologies

Omni-Lok™ is a trademark of Omnifit Ltd.

Silicone Select™ is a trademark of Bio-Chem Valve Inc.

www.beionmed.com (提供中文下载服务)
021-51096593 Email: sales@beionmed.com

Rev. 0205