



Connection Kits



Everything you need to make hardwall tubing connections

- Inert connectors and fittings
- Omnifit Gripper fittings rated to 1000psi (68 bar)
- Kits for 1/16" (1.6mm) and 1/8" (3.2mm) outside diameter tubing

Tube Connector Kits

These kits contain parts to make 10 Gripper fittings and come complete with 5 coupling bodies to enable simple straight-through connection of two Gripper fittings. Available in 1/16" and 1/8" OD tubing versions and with either polypropylene or Tefzel® couplings.

Starter Kits

Omnifit's starter kits provide a range of fittings and adaptors selected to give a variety of tube connection options for the user. Each kit contains enough Omnifit Gridders, fittings and PTFE tubing to make 20 fittings. Also included are 10 x 2-way couplings, 1 male and 1 female Luer fitting and 1 x 3-way 'T'-connector. Starter kits are available in 1/16" and 1/8" OD tubing versions.

Laboratory Kit

A comprehensive kit of parts to enable the user to make a range of tubing connections with versions accommodating tubing up to 11mm O.D. Three different valve and connector types provide a wide range of flow splitting and distribution options.

Gripper Trial Pack

The Gripper Trial Pack contains 2 gripper fitting nuts, 2 PTFE gridders, a Tefzel® coupling and 1/8" O.D. PTFE tubing.

Ordering Information

Tube Connector Kits

1/16" and 1/8" Polypropylene
Tube Connector Kits
Part Numbers 2003 and 2013

Contents:

- 10 x polypropylene tube end fittings (2110, 2210)*
- 5 x 2-way polypropylene couplings (2301)*
- 10 x PTFE Grippers (2310, 2312)



2003



2013

1/16" and 1/8" Tefzel®
Tube Connector Kits
Part Numbers 2004 and 2014

Contents:

- 10 x Tefzel® tube end fittings (2130, 2230)
- 5 x 2-way Tefzel® couplings (2302)*
- 10 x PTFE Grippers (2310, 2312)



2004



2014

Starter Kits

1/16" Starter Kit
Part Number 2002

Contents:

- 20 x tube end fittings (2110)*
- 10 x 2-way couplings (2301)
- 1 x female Luer fitting (2502)
- 1 x 3-way 'T'-connector (2401)
- 20 x Grippers (2310)*
- 1 x male Luer fitting (2502)
- 2 x plugs (2320)
- 20m x 1/16" OD PTFE tubing (008T16-080-20)



2002

1/8" Starter Kit
Part Number 2012

Contents:

- 20 x tube end fittings (2210)*
- 10 x 2-way couplings (2301)
- 1 x female Luer fitting (2502)
- 1 x 3-way 'T'-connector (1010)
- 20 x Grippers (2312)*
- 1 x male Luer fitting (2502)
- 2 x plugs (2320)
- 20m x 1/8" OD PTFE tubing (008T32-150-20)*



2012

Laboratory Kit

Part Number 2015

Contents:

- 2 x 2-way connector (1001)
- 1 x universal connector (1008)
- 1 x 2-way 'T'-valve (1109)
- 10 x 1/16" tube end fittings (2110)
- 10 x 1/16" Grippers (2310)
- 1 x 3-way 'T'-connector (1010)
- 1 x distribution valve (1118)
- 2 x 2-way coupling (2301)
- 20m x 1/16" OD/0.8mm ID PTFE tubing (008T16-080-20)



2015

Gripper Trial Pack

1/8" Gripper Trial Pack
Part Number 2011

Contents:

- 2 x tube end fittings
- 2 x Grippers
- 1 x 2-way coupling
- 2m x 1/8" OD/1.5mm ID PTFE tubing



2011

* Colors or quantity may differ from that sold as stand alone product.

Trademarks:

Tefzel® is a registered trademark of E.I. du Pont de Nemours and Company

Bio-Chem Valve / Omnifit, 2 College Park, Coldhams Lane, Cambridge, CB1 3HD UK
t: +44 (0) 1223 416642 f: +44 (0) 1223 416787 e: sales@omnifit.com w: omnifit.com

Rev. 0205

Bio-Chem Valve / Omnifit, 85 Fulton Street, Boonton, NJ 07005 USA
t: 973 263 3001 f: 973 263 2880 e: info@bio-chemvalve.com w: bio-chemvalve.com